





PRODUCT DATA PACK (PDP)	PDP REF: PDP-ANIGENE WIPE V01	
NAME: ANIGENE Professional Surface Disinfectant Wipes CODE: B3015F		

DOCUMENT ISSUE			
ROLE	NAME	JOB TITLE	SIGN & DATE
AUTHOR	Matthew Hunt	Microbiologist	<u>Matthew Hunt</u> <small>Matthew Hunt (Jun 26, 2024 11:29 GMT+1)</small>
REVIEWER	Dr Chris Plummer	Head of Technology	<u>C. Plummer</u> <small>C. Plummer (Jun 26, 2024 11:31 GMT+1)</small>
QUALITY APPROVER	Megan Bowles	Quality Specialist	<u>Megan Bowles</u> <small>Megan Bowles (Jun 26, 2024 11:57 GMT+1)</small>
ISSUED BY	Megan Bowles	Quality Specialist	<u>Megan Bowles</u> <small>Megan Bowles (Jun 26, 2024 11:57 GMT+1)</small>
ISSUE DATE	26 JUN 2024	REVIEW DATE	26 JUN 2026

PRODUCT DATA PACK (PDP)	PDP REF: PDP-ANIGENE WIPE V01	
NAME: ANIGENE Professional Surface Disinfectant Wipes CODE: B3015F		

PRODUCT INFORMATION	
PRODUCT NAME	ANIGENE Professional Surface Disinfectant Wipes
PRODUCT CODE	SKU100A
FORMULATION VARIANTS	B3015F - Citrus
SDS REFERENCE	SDS-B3015F ANIGENE Professional Surface Disinfectant Wipes_English_v1.0.0
<b>NOTE: ACTIVE INGREDIENTS AND SAFETY INFORMATION DETAILED IN SDS</b>	
DIRECTIONS FOR USE	<p>Suitable for use on hard, non-porous surfaces. Lift seal and remove one wipe from pack. Place palm of hand on top of the wipe. Apply the wipe to the surface in an 'S' pattern, ensuring full coverage. Each wipe is suitable for disinfecting 1m<sup>2</sup>. Avoid wiping the same area twice and use a new wipe for each new surface, to avoid cross-contamination. Leave surface to dry naturally, no need to rinse. Test an inconspicuous area before use if unsure. Reseal pack after use to avoid wipes drying out. Always rinse food / water bowls and toys after disinfection.</p> <p>Dispose of used wipes in general waste or clinical waste as appropriate. Wipes are biodegradable within 3 months.</p>
APPLICATION AREAS	PT 3 – Veterinary Hygiene
STORAGE CONDITIONS	<p>Keep out of reach of children. If medical advice is needed, have product container or label at hand. Avoid release to the environment. Dispose of contents/container in accordance with national regulations. Incompatibility with other materials: Oxidising Substances, Anionics, Hypochlorite, Aldehydes.</p> 
SHELF-LIFE	<p>2-Years</p> <p>Established following criteria for accelerated storage test in section 3.4.1.1 of ECHA guidance on the BPR: Volume I. Part A. Chapter III: Requirements for Biocidal Products. Version 1.1 November 2014</p>


PRODUCT SPECIFICATION	
APPEARANCE	Colourless to pale yellow impregnated wipes
COLOUR	Colourless to pale yellow
ODOUR	Characteristic
pH	9.00 – 10.00
SPECIFIC GRAVITY	0.990 – 1.020
VISCOSITY	N/A

<b>PRODUCT DATA PACK (PDP)</b>	<b>PDP REF: PDP-ANIGENE WIPE V01</b>	
<b>NAME: ANIGENE Professional Surface Disinfectant Wipes CODE: B3015F</b>		

BACTERICIDAL EFFICACY					
TEST REF	TEST ORGANISMS	TEST CONDITIONS	TEST RESULT	REPORT REF	COMMENTS
EN1656	<i>Enterococcus hirae</i>	1 min, low soil, 20°C	>5 log reduction	TR-910	N/A
EN16615	<i>Enterococcus hirae</i>	1 min, low soil, 20°C	>5 log reduction	TR-913	Modified stainless steel
EN1656	<i>Pseudomonas aeruginosa</i>	1 min, low soil, 20°C	>5 log reduction	TR-910	N/A
EN16615	<i>Pseudomonas aeruginosa</i>	1 min, low soil, 20°C	>5 log reduction	TR-913	Modified stainless steel
EN1656	<i>Staphylococcus aureus</i>	1 min, low soil, 20°C	>5 log reduction	TR-910	N/A
EN16615	<i>Staphylococcus aureus</i>	1 min, low soil, 20°C	>5 log reduction	TR-913	Modified stainless steel
EN1656	<i>Proteus hauseri</i>	1 min, low soil, 20°C	>5 log reduction	TR-911	N/A
EN16615	<i>Proteus hauseri</i>	1 min, low soil, 20°C	>5 log reduction	TR-913	Modified stainless steel
EN1656	<i>Escherichia coli</i>	1 min, low soil, 20°C	>5 log reduction	TR-911	N/A
EN1656	<i>Brucella intermedia</i>	1 min, low soil, 20°C	>5 log reduction	TR-926	N/A
EN1656	<i>Bordetella bronchiseptica</i>	1 min, low soil, 20°C	>5 log reduction	TR-926	N/A
EN1656	<i>Streptococcus equi</i>	1 min, low soil, 20°C	>5 log reduction	TR-926	N/A
EN1656	<i>Salmonella typhimurium</i>	1 min, low soil, 20°C	>5 log reduction	TR-927	N/A
EN1656	MRSE ( <i>Staphylococcus epidermidis</i> )	1 min, low soil, 20°C	>5 log reduction	TR-927	N/A
EN1656	MRSP ( <i>Staphylococcus pseudointermedius</i> )	1 min, low soil, 20°C	>5 log reduction	TR-927	N/A

YEASTICIDAL/FUNGICIDAL EFFICACY					
TEST REF	TEST ORGANISMS	TEST CONDITIONS	TEST RESULT	REPORT REF	COMMENTS
EN1657	<i>Candida albicans</i>	1 min, low soil, 20°C	>4 log reduction	TR-912	N/A
EN16615	<i>Candida albicans</i>	1 min, low soil, 20°C	>4 log reduction	TR-914	Modified stainless steel
EN1657	<i>Aspergillus brasiliensis</i>	30 min, Low soil, 20°C	>4 log reduction	TR-921	N/A
EN16438	<i>Aspergillus brasiliensis</i>	30 min, Low soil, 20°C	>3 log reduction	TR-922	N/A

VIRUCIDAL EFFICACY					
TEST REF	TEST ORGANISMS	TEST CONDITIONS	TEST RESULT	REPORT REF	COMMENTS
EN14675	Vaccinia virus	5 min, low soil, 20°C	>4 log reduction	TR-917	N/A
EN17122	Feline coronavirus	5 min, low soil, 20°C	>4 log reduction	TR-918	N/A

<b>PRODUCT DATA PACK (PDP)</b>	<b>PDP REF:</b> <b>PDP-ANIGENE WIPE</b> <b>V01</b>	
<b>NAME: ANIGENE Professional Surface Disinfectant Wipes</b> <b>CODE: B3015F</b>		

ADDITIONAL DATA		
TEST METHOD	RESULT	TEST REPORT REF
Preservative Efficacy Testing	B3015 complies with the test for the Efficacy of Antimicrobial Preservatives (European Pharmacopoeia 11 <sup>th</sup> Edition)	TRC-338

CHANGE HISTORY			
VERSION	CHANGE REFERENCE	CHANGE TYPE	REASON FOR CHANGE
01	DCC24.078	NEW	New PDP for new product











# PDP-B3015F-ANIGENE v01

Final Audit Report

2024-06-26

Created:	2024-06-26
By:	Megan Bowles (mbowles@byotrol.com)
Status:	Signed
Transaction ID:	CBJCHBCAABAAB2_ypSn86DLXGmZyWgQiV-VF7zLSGT7M

## "PDP-B3015F-ANIGENE v01" History

-  Document created by Megan Bowles (mbowles@byotrol.com)  
2024-06-26 - 9:08:36 AM GMT
-  Document emailed to mhunt@byotrol.com for signature  
2024-06-26 - 9:08:40 AM GMT
-  Email viewed by mhunt@byotrol.com  
2024-06-26 - 10:28:52 AM GMT
-  Signer mhunt@byotrol.com entered name at signing as Matthew Hunt  
2024-06-26 - 10:29:29 AM GMT
-  Document e-signed by Matthew Hunt (mhunt@byotrol.com)  
Signature Date: 2024-06-26 - 10:29:31 AM GMT - Time Source: server
-  Document emailed to cplummer@byotrol.com for signature  
2024-06-26 - 10:29:32 AM GMT
-  Email viewed by cplummer@byotrol.com  
2024-06-26 - 10:29:57 AM GMT
-  Signer cplummer@byotrol.com entered name at signing as C. Plummer  
2024-06-26 - 10:31:24 AM GMT
-  Document e-signed by C. Plummer (cplummer@byotrol.com)  
Signature Date: 2024-06-26 - 10:31:26 AM GMT - Time Source: server
-  Document emailed to Megan Bowles (mbowles@byotrol.com) for approval  
2024-06-26 - 10:31:28 AM GMT





Document approved by Megan Bowles (mbowles@byotrol.com)

Approval obtained using URL retrieved through the Adobe Acrobat Sign API

Approval Date: 2024-06-26 - 10:57:49 AM GMT - Time Source: server



Agreement completed.

2024-06-26 - 10:57:49 AM GMT



**Adobe Acrobat Sign**