



A world of potential

VET Capnography Application Suggestions



EDAN Instruments, Inc.

Luca.zuo@edan.com

www.edan.com



2024.04.11



CO2 Configuration

Surgery setting

Click CO2 reading area to enter **CO2 Setup** >> **Other Setups**



For surgery, normally a mixture of pure oxygen and a little AG are being used. To ensure the accuracy of the measurement, it's required to set the O2 compensation to **100%**, Anesthesia Agent compensation to **around 1.5%**

- The rest parameters are unchanged as default.
- For very small animals, in case it couldn't give an ideal EtCO2 reading, we could try to decrease the pump rate.



General issues

Three major issues in practices

**EtCO2
no
reading**

**EtCO2
low
reading**

**FiCO2
high
reading**

● Troubleshooting

Is there an alarm message?



Airway is occluded.

- Check the status of the watertrap and line.
- try to replace them
- Check internal connectors/tubes

Two possible causes

1. If the above **CO2 Occlude** error is not solved for a long time, then the system will show **CO2 Sensor Faulty** and stop taking in air as a protection action.
2. The G2 module is indeed faulty. Replacement is needed.

● Troubleshooting

EtCO2 no reading—without alarm message

Try to breathe into the sampling line, see if it shows any reading and trace

● Troubleshooting

EtCO2 no reading—without alarm message

If it WORKS with human breathe

1. Verify If the animal is breathing

For some cases the patient is barely breathing, so it could result in no reading. In some cases, changing the position of the patient may help

2. Verify the intubation is correct


watch out for oesophageal intubation

3. Obstruction

ET tube is obstructed, try to verify it, install it again or move it a little bit

4. Air leakage

- Make sure the size of ET tube is suitable, the patient must not breathe around the tube



- Make sure all connections are firm

● Troubleshooting

EtCO2 no reading—without alarm message

If it **DOES NOT WORK** with human breathe

1. Replace the watertrap and line

It being leaking or blocked could be the cause, try to replace them with new ones.
Watertrap and sample line are disposable and must be replaced for each patient.

2. Make sure the watertrap is firmly installed

Watch out for possible loose contact which may cause leakage

3. Make sure the line is firmly installed on the watertrap

Watch out for possible loose contact which may cause leakage

4. Check the bracket O ring


The rubber O ring could be worn out with time and cause leakage. Needs to be changed if you notice any damage





Troubleshooting

EtCO2 LOW reading—without alarm message




● Troubleshooting

EtCO2 LOW reading—without alarm message

Other causes

1. Dead space


- a. Use ET tube with side-port connector



- b. Use low dead-space ET tube adaptor with luer sampling port



- c. Make sure that the ET tube is not excessively long



- d. For small patient, do not use L tube (right angle)
e. Sample close to ET tube as possible

2. Flow rate of oxygen high

It could interfere and dilute the CO2 concentration.
Try to decrease it

3. Module needs calibration

Send it back for service





Troubleshooting

EtCO2 LOW reading—without alarm message

Try to breathe into the sampling line, see if it shows correct reading, around 38mmhg

Troubleshooting


FiCO₂ high reading—without alarm message



● Troubleshooting

Any abnormal reading or trace

Try zero calibration



Note:

Calibration must be done in a ventilated environment.


It's very important to make sure the space is well-ventilated during the surgery. Because auto zero calibration is being done from time to time.

● In case problem cannot be solved

In case the problem cannot be solved with the suggestions mentioned above,

Please try to take a **VIDEO** in case of any problem.

A video including the CO2 trace, the CO2 reading, the tubes, the connections and the ET tube would be very helpful for further troubleshooting.






Standby mode

When a surgery is finished, remember to put the monitor to **STANDBY MODE** in terms of EtCO₂ function or just switch off the monitor.

Click CO₂ reading area to enter **CO₂ Setup**

>> **Work Mode**

>> select **Standby**



In standby mode, the EtCO₂ module stops working. The pump rests and it stops sampling.

Benefit:

- It avoids dust and moisture from entering the system.
- It helps to extend the life-span of the module

Resume measurement:

To start the CO₂ measurement, select **Measure** for work mode



Quick summary

Always consider replacing the watertrap and line

Try to reduce dead space-particularly for small patients

Try to test with human

Try auto zero

Compensation settings must be done

Put the monitor to Standby mode

Try to take a video in case of any problems. The video should include the trace and reading, the tubes and connections, ET tube



A world of potential

THANK YOU

Edan Instruments, Inc.