

Vet Direct

A Great Deal More For Less!!

Ajex 140H

- For Stress Free Radiology -



FEATURES:

- ◆ Easy to Use – Error Proof Digital Control Panel
- ◆ Superb High Frequency Power
- ◆ Up to 100kV
- ◆ 2.4kW
- ◆ Advanced 60kHz High Frequency Generator
- ◆ Minimal Exposure Times
- ◆ Compact
- ◆ Long Life Toshiba X-ray Tube
- ◆ High Quality X-ray with Superb Resolution
- ◆ State-of-the-Art Design

HIGH FREQUENCY GENERATOR

Ajex 140H offers your practice a powerful high frequency generator for stunning high quality x-rays with reduced exposure time.

SPECIFICATIONS

High Frequency Generator:	60KHz
Focal Spot:	1.2mm
Max. Output	2.4kW
Rating:	50~100kV 40mA~15mA
kV:	Continuous 50 to 100kV (2kV Steps)
mA:	5mA Steps
Exposure Time:	0.8mAs - 120mAs (0.02 - 3.0 seconds)
X-Ray Control:	Hand Switch
Tube:	Toshiba 1.2mm / D124
Exposure Control:	All parameters are controlled by microprocessor
Dimensions:	238 (W) x 340 (D) x 194 (H) mm
Weight:	14.5Kg

MICROPROCESSOR CONTROLLED PANEL

The Control Panel on the Ajex140H is the "brains" of the system: all exposures are controlled by a powerful microprocessor.

The controls are simple, intuitive and adaptable. The control unit boasts large clear "Soft Touch" keys. There are 6 programmable keys for your most commonly used vectors.

Price

Code	Description	Price
AJEX140	Ajex 140H X-ray	£5600

Finance Estimate

Ajex140H	36 Months		60 Months	
Amount	Monthly	Daily	Monthly	Daily
£5600	£180	£8.27	£115	£5.30

The rentals above are subject to VAT, which is reclaimed by the practice in normal returns. Daily cost is based on 5-day week. (No allowance deducted for tax relief). Rentals are fully allowed as a business expense for tax purposes and claimed annually. At the end of the agreement the use of the equipment is available for a secondary fixed annual rental of £25+VAT. Full written details available on request. E&OE

Prices Correct 01/10/04 E&OE. All prices exclude VAT