

Weck Horizon Ligation Clips

Dimensions – Width, Height and Closed Length

Clip Description Sorted by clip type.	Width		Height		Closed Clip Length	
	Inches	Millimetres	Inches	Millimetres	Inches	Millimetres
White Micro Horizon	0.075	1.91	0.075	1.91	0.088	2.24
Small Red Horizon	0.082	2.08	0.113	2.87	0.145	3.68
Small Yellow Horizon	0.078	1.98	0.113	2.87	0.143	3.63
Medium Horizon	0.119	3.02	0.184	4.68	0.232	5.89
Med-Large Horizon	0.226	5.73	0.284	7.21	0.353	8.98
Large Horizon	0.307	7.80	0.396	10.07	0.488	12.38

