

The Art of Kennelling



It's the **little** things we do
that make a
BIG difference



Kennelling is a vital part of any modern-day practice. Patient wellbeing during the recovery process is more important than ever.

At Vet Direct we are proud to introduce the Stellvet range. A brand synonymous with innovative veterinary driven solutions where patient comfort and ease of use are the primary goals and offer outstanding quality at an affordable price.

This all new range of kennels, accessories and walk in kennels has been developed from first principles. Development from client feedback and market analysis has produced features that overcome problems with conventional designs and bring distinct benefits to users and patients alike.

With the modular kennels, sizes are comparable with other industry leading brands and the system is compatible with many existing kennel ranges.

Walk in kennel doors feature an all new glass sealing design that take hygiene and serviceability to a new level.


Attention to detail is everywhere and more than anything the Stellvet range really does reflect our motto 'It's the little things we do that make a big difference'

We really believe that you will see, feel and experience a noticeable difference with the new Stellvet range.

THE STELLVET DIFFERENCE

Designing an all new kennel range is no small task, especially one that incorporates many new features and addresses issues that have plagued other manufacturers for many years.

At Vet Direct we do not compromise and believe what has been created offers the finest solution available to the discerning veterinarian today. Don't just take our word for it, experience the difference that is Stellvet for yourself.



Soft Close Latches



Integral Sealing Strip



Quick Release Gate




KENNEL SIZES

All Stellvet kennels are manufactured from 304 grade stainless steel. Material thickness varies depending on size of the kennel selected. All kennels are compatible with market leading brands and can be integrated, using Stellvet connectors, with existing banks from other manufacturers.

All sizes are shown in the table below and are given as Width x Height x Depth in millimetres (metric) in the first column and inches (imperial) in the second column.

W x H x D (mm)	W x H x D (inches)	Single Door / Double Door
457 x 457 x 558	18 x 18 x 22	Single Door
610 x 457 x 558	24 x 18 x 22	Single Door
610 x 457 x 717	24 x 18 x 28	Single Door
610 x 610 x 717	24 x 24 x 28	Single Door
610 x 762 x 717	24 x 30 x 28	Single Door
610 x 914 x 717	24 x 36 x 28	Single Door
762 x 610 x 717	30 x 24 x 28	Single Door
914 x 610 x 717	36 x 24 x 28	Single Door
762 x 762 x 717	30 x 30 x 28	Single Door
762 x 914 x 717	30 x 36 x 28	Single Door
914 x 762 x 717	36 x 30 x 28	Single Door
914 x 914 x 717	36 x 36 x 28	Single Door
1220 x 762 x 717	48 x 30 x 28	Double Door
1220 x 914 x 717	48 x 36 x 28	Double Door



Sizes given are the maximum external dimensions of the kennels. Single door kennels all feature a reversible wire mesh door that can be hinged from the left of the right-hand side easily using the mounting hardware provided. Double door kennels feature no latch bar to provide uninterrupted patient access.

ACCESSORIES



A wide range of accessories are available to individually tailor your Stellvet kennels to your exact requirements. These include:

Fixed Platforms	Mobile Platforms	Side Panels	Back Panels
Top Panels	High Level Platforms	High Level Cupboards	Low Level Cupboards
Glass Doors	Heating	Electrical Rails	Service Rails
Drawer Units	Clipboards	Ventilation & Extraction	Integrated Desks
Lighting	IV Holders	Infusion Pump Holders	Medicine Holders

COLOURED LATCHES

A wide variety of coloured latches are available to suit any interior design. We even have two-tone options. Colours include:

Black **White** **Red** **Yellow** **Green** **Blue** **Light Grey** **Dark Grey**



WALK IN KENNELS

Accommodating larger dogs very often requires specialist walk in kennels. As well as a modular kennel range, Stellvet offer a wide selection of walk in kennel doors, complete pen systems and catteries.



Doors, side panels, dividing walls and ventilation grills are available in a wide range of styles including glass, bars, mesh, trespa or combination units.

Products can be custom made to any bespoke size but as a standard the doors are produced to dimensions of 800mm width and 2000mm height.



THE COMPLETE SERVICE

At Vet Direct we believe in providing our customers with the finest service. The Stellvet range isn't just manufactured, it is designed, planned, installed and supported.

Pre-installation checks are carried out by experienced surveyors and all of the installation works undertaken by a dedicated team of site fitters, all time served craftsmen.

In addition, should your existing kennel systems require modification or repair, we will be happy to assist.

If you have any requirements for new equipment, a major new project or simply a basic repair, we will treat any enquiry with the utmost of importance and be only too happy to help.

We have one goal in mind – total customer satisfaction.





Contact details

Freephone: 0800 068 3300

Freefax: 0800 068 3311

Int Tel: +44 (0) 191 264 2705

Int Fax: +44 (0) 191 268 3311

Email: info@vet-direct.com

Twitter: @vetdirectuk

Website: www.vet-direct.com