

Edan M3B

Capnograph with Pulse Oximeter



Image Shown is Edan M3B with Sidestream Module Attached

Edan M3B

With SpO2 (Nellcor optional) and Respiration CO2 highly advanced monitoring technologies, EDAN M3B ensures effective capnography monitoring for intubated and non-intubated patients for continuous long-term monitoring.

Available in 2 Configurations:

- SideStream EtCO2 & SpO2
- MainStream EtCO2 & SpO2



Features:

5.7" inch Colour TFT Screen

Lightweight, Portable Design and User-Friendly Interface for Easy Operation

Backlight Control and Standby Mode for Power Saving

Respironics Sidestream OR Mainstream CO2 Measurement

Display Numeric and Waveform Information Simultaneously

Powerful Storage Capacity

Built-in Lithium-ion Battery for long continuous work (Up to 10 Hours)

Suitable for Adult, Paediatric and Neonatal Patients

Trend Table Review and Record

Trend Graph Review and Record

USB Data Storage and Review

Size and Weight:

173 (L) x 241 (H) x 189 (D) mm

Battery:

Battery Lithium 14.8V 4,400mAh

Working period 600 min

Recharge period 300min

Mounting Solutions:

- Wall mount
- Rolling Stand
- IV Pole Clamp

