

## Mag-Tek Binocular Microscope XSB2010



This is a great quality microscope at an economy price.

Features:

- Siedentopf type binocular head, inclined at 45 degrees
  - Interpupillary distance: 55-75mm
  - Eyepieces: WF 10X
  - Objective: Quadruple mechanism
- Achromatic Objective Lenses: x4, x10, x40x 100 (oil)
- Stage: Double layer mechanical stage with removable slide holder
  - Focusing: Coaxial coarse and fine focussing
- Condenser: Abbe NA 1.25 condenser with iris diaphragm
  - Illumination: Halogen 6V 20W

**Vet Direct Price £295**