

Mag-Tek Binocular Microscope XSB2030 with Trinocular



Price £555.90

Biological Microscope with Finitiy Seidentofp Binocular Head 30 degree inclined, 360 degree rotation.

- * Trinocular Head for Fitting Camera
- * Magnification: 40X~1000X
- * Eyepiece: WF10X, FOV=20mm
- * Interpupillary distance: 50~75mm, diopter adjustment
- * Achromatic Objectives(4pcs): 4x, 10x, 40x(s,oil), 100x(s,oil)
- * Stage: Mechanical stage with removable slide holder 142x140mm; Travel area 76(X)x50(Y)mm with a right-hand stage handle; 0.1mm provide; Double slide holder;
- * Condenser: Abbe type N.A.=1.25 with Iris Diaphragm.
- * Focusing: Coaxial coarse/fine focusing knobs.
- * Nosepiece: Inward Quadruple click-stop - Revolving mechanism with multiple ball bearing
- * Illumination: Built-in Koehler illuminator system, halogen 6V-20W, Adjust brightness, Field Diaphragm

Mag-Tek Microscope Camera Eyepiece 3.1MP Full Kit

USB Microscope Camera (with software) 3.1M pixels Camera Eyepiece for PC with USB & software



Price £245.00