

Pulse Oximeter H100

Edan



The VE-H100B Veterinary Pulse Oximeter

Featuring:

- LCD display
- RTC (Real Time Clock) display
- Backlight switch and automatic power-off function for power saving
- SpO2 and Pulse rate measurement
- Display number and waveform of SpO2 simultaneously
- Long battery life - up to 48 hours on 4 "AA" or Ni-H rechargeable batteries
- Powerful data storage capacity
- Audible and visible alarm capability - Pitch Tone

Provides a simple universal tool for veterinary use – from spot checks to continuous monitoring.

Bright and clear LCD Screen offers you SpO2 and Pulse Rate measurement, the information can be displayed as numeric and waveform and simultaneously.

The latest 10 minutes trend graph and table of SpO2 and Pulse Rate can be reviewed on-screen on the fly.

Convenient size that will work up to 48 hours on AA batteries – Internally it can hold 100 IDs and 300 hours data storage capability for continuous monitoring.

Standard Accessories:

Main Unit
EDAN Veterinary Pulse Oximetry Sensor
AA Size Batteries
Carrying Case
User Manual VE-H100B
Quick Reference Card

Optional Accessories:

Desk Top Charger
Protective Cover
Data pack*

*Data can be transferred to PC for storage, review and printing - You will need data transfer lead and Edan software to allow data to be accessed and downloaded.