


DeltaWave® 400

Anaesthetic System Complete



The DeltaWave® 400 Complete Comprises Of:

1. **DeltaWave® 400**

Non-Rebreathing Anaesthesia System is designed for safe, simple, and controllable inhalation. DeltaWave® adapts to the needs of your surgery. Oxygen rotameter featuring a glass flowmeter from 0-10L/min Accepts Sevoflurane, Isoflurane or Halothane cagemount vaporisers. Small foot print size.

2. **New Tec3 Vaporiser**

The Tec Type3 vaporiser has a simple, yet durable design, excellent performance and high degrees of safety and reliability.

As these are NEW units they do not have to be serviced or calibrated for 3 years from purchase. Meaning a saving of two complete service and calibrations over buying a reconditioned unit, saving on average £300+ in engineer costs in addition to the hassle and time.

3. **VetO2 Plus Oxygen Concentrator**

10L Single Flow Filters ambient air and produces 0 to 10L single flow of oxygen per minute Robust body, easy to clean Compatible with most anaesthesia systems
Features: flow control and timing display, alarm

The DeltaWave® Standard Floor Model is our most popular machine.

- Contemporary Design
- Small Foot Print
- Easy Roll Stand with braked wheels
- Oxygen Flush
- Anaesthetic Circuit Holder
- Storage Shelf
- Easy to clean