

DeltaWave 200

Trolley Anaesthetic System



The DeltaWave® Non-Rebreathing Anaesthesia System designed for safe, simple, and controllable inhalation. DeltaWave® adapts to the needs of your surgery.

- Non-rebreathing anaesthetic delivery system
- Oxygen rotameter featuring a glass flowmeter from 0-10L/min
- Accepts Sevoflurane or Isoflurane cagemount vaporisers

The DeltaWave® Standard Floor Model is our most popular machine.

- Contemporary Design
- Small Foot Print
- Easy Roll Stand with braked wheels
- Accepts both E & F Oxygen Cylinders
- Holds Two Cylinders
- Oxygen Flush
- Anaesthetic Circuit Holder
- Monitor Shelf
- Easy to clean

Optional Extras

Vaporiser (Sevoflurane, Isoflurane)

DeltaWave® Dual Rotameter 0-1L and 1-10L Cascade for Small Patients*

DeltaWave® Dual Rotameter Oxygen and Nitrous

Mayo Tray/Monitor Shelf

IV Pole


Low Oxygen Alarm

Circle Converter Adapter

Compatibility block for cagemount vaporisers

Oxygen Regulator

Note: Price excludes circuit, vaporiser and fluosorber depicted on image.



Also available in tabletop and wall mount configurations