

Mini Centrifuge



Our unique design allows for the interchangeable rotors on our centrifuge to be converted from centrifuge tube to strip tube in just seconds.

The Mini Centrifuge is supplied with two rotors, one for microtubes and one for 0.2ml PCR strip tubes. The 6 x 1.5ml rotor is designed to centrifuge up to six individual 1.5ml plastic micro centrifuge tubes. It will also accept 0.5ml and 0.4ml tubes with the adapters supplied with the unit. Adapters are available separately for 0.2ml tubes. The strip tube rotor is designed to centrifuge two 8 x 0.2ml strip tubes. It will also accept up to 16 individual 0.2ml tubes. Both rotors are designed for applications requiring relatively low g-forces, such as micro filtration; cell separation and quick spin downs from the walls of tubes. Please read this manual thoroughly before attempting to operate the centrifuge.

General Features:

- Interchangeable rotors Accommodate 1.5 and 0.5 Tubes and 0.2 ml Strip Tubes
- Durable hinge
- Simple Latching device for easy access
- Safe - Activate the unit without having to open the Lid
- ISO 9001 Registered manufacturer

Includes:

One centrifuge tube rotor for 1.5 ml tubes

Adapters for 0.5 ml tubes

One strip tube rotor

One locking key

Specifications:

Speed: 6000 RPM

Dimensions: W: 126mm, H: 100mm, D: 138mm