

# Op / Prep Table

OPTVD01



## Operating / Prep Table Hydraulic

- Stainless Steel Top operation / prep table with reinforced construction that eliminates wobbling and tipping
- Top is dished with drainage hole and removable and easily washable collection basin
- Top has plenty of cleats for tie-down
- Includes 110cm height adjustable IV pole
- Easy clean top with drain and basin
- Table top size 127cm x 66cm
- Easy foot operated Hydraulic pump height adjusts from 56cm-101cm high
- Adding wheels to this system will add around 10cm to the height bandwidth

Code: OPTVD01