

# DenTek D5

Piezo Scaler with LED



LED with light effect



Without light effect

## D5 LED

DenTek-Excel D5 Piezo Scaler

### Features:

- Titanium alloy handpiece with auto tuning piezo technology
- LED Light in handpiece – to illuminate scaling area
- Autoclavable detachable handpiece to effectively avoid cross-infection
- Programmable Buttons to Store: General Scaling, Perio and Endo Functions
- Cavitations produced by the scaling tip ensures easy scaling and also kills germs
- Automatic frequency tracking to ensures your system is optimised
- Small and light design for easy operation: 188mm x 159mm x 81mm
- Easy to control water spray to ensure tip remains cool

Superb quality and value

The system is compatible with SATELEC's handpieces & tips

12-Month Warranty

Includes:

Five scaling tips  
Detachable titanium alloy handpiece  
Spare LED light ring  
Foot pedal  
Power supply  
Water connection  
Water supply  
Scaling tip wrench

Automatic degenerative feedback, constant power output, and more comfortable scaling; Detachable handpiece with LED light; Functions: Scaling, Perio, Endo

### **Technical Parameters**

Main unit input	24V-50Hz/60Hz 1.3A
Output power	3W-20W
Frequency	28kHz $\pm$ 3kHz
Water pressure	0.1bar ~ 5bar(0.01MPa ~ 0.5MPa)
Main unit weight	1kg
Dimensions	188mm×159mm×81mm