

Practical Trolley 49 x 37 x 77H cm



Practical Trolley 49 x 37 x 77H cm (Delivered in Kit Form)

A practical trolley for ECG, Electrosurgical unit or any other device, with 4 nylon castors and 2 lateral handle for easy moving.

Technical Specifications

- Revolving castors Ø 50 mm
 - Size: 49 x 37 x h 77 cm
 - Shelf size: 40 x 30 cm
 - Colour: RAL 9002
 - Delivered in kit form

(NOTE: All kit on trolley image is not included)

Vet Direct unit 9-10 Gateway west Kingfisher Boulevard Newburn Riverside Newcastle Upon Tyne
NE15 8NZ

Tel: 0800 0683300 Fax: 0800 0683311 Int: +44 (0)191 2683300
Email: info@vet-direct.com Web: www.vet-direct.com

PT27438 £125

Vet Direct unit 9-10 Gateway west Kingfisher Boulevard Newburn Riverside Newcastle Upon Tyne
NE15 8NZ

Tel: 0800 0683300 Fax: 0800 0683311 Int: +44 (0)191 2683300
Email: info@vet-direct.com Web: www.vet-direct.com