

Vet Direct

...on the pulse

Cozy™ Warm

Patient Warming System by Darvall



Cozy™ Warm - Air Warming System New Unit comes with 24 Month Guarantee

Fast warm-up time just 75 seconds

Features:

1. Advanced interface with 4 one-touch temperature selections (Ambient, 32°C, 38°C, 43°C) displayed on large digital readout to meet all therapy needs;
2. Hose-end temperature sensing ensures accurate air temperature delivery
3. Small size frees up more usable workspace
4. Light weight makes transport and set-up a breeze
5. Low noise. The quiet Warming Unit can be placed conveniently under the operating table, head or foot end.

Safety first:

1. Built-in hour meter makes it easy to monitor usage for preventative maintenance
2. Over/under temperature and disconnect audible and visual alarms ensure safety
3. High efficiency 0.2 micron air filter reduces airborne particulate and ensures clean air delivery

Will fit with "Bair Hugger" style blankets and accessories



Read on for comparisons "v" bair hugger

Comparison of Vet Force with Bair Hugger

Attributes	Cozy Warm	Bair Hugger
Model	#9612	505; 750
Types of Warming Unit	Only one type.	Two types available – one for large blankets only and one for smaller blankets only.
Dimension	28cm*22cm*22cm	33cm*28cm*25cm
Weight	4.5kg	6.2kg or 7.4kg
Air Flow	High Level: 48 CFM Low Level: 35 CFM	28-30 CFM
Heater	PTC	Resistance Wire
Heater Power	1200w	850w; 1400w
Time to Reach the Setting Temperature	75s	More than 3 minutes
Filter	0.2 µm	0.2 µm
Filter Life	500 hours	500 hours
Noise Level	High Level<55db Low Level<50db	55db
Temperature Sensor Number	2	2
Waterproof Grade	IP21	IP21
Displayer	LED	Digital Tube Display
Price	Competitive	Expensive
Range of blankets	Very wide range	Very wide range
Blanket compatible with other brand warming unit	Yes	No
Warming Unit Warranty	2 years	

