

NEW Laryngoscope Fits to your standard handle



The NEW and completely unique Opticlar Veterinary Laryngoscope Range provides the Veterinary Professional with superb quality LED laryngoscope instruments that fits to their existing Otoscope / Ophthalmoscope handles (For Opticlar, WA, Heine or Keeler) Offering you the following benefits:

- Superior and brighter illumination over standard laryngoscopes
- Magnifier Option makes this a very useful all round tool*
- No more changing bulbs (Or hot bulbs)
- No more discharged or corroded batteries
- 10 Year LED Warranty
- 10 Year Physical Product Warranty
- Uses existing Oto/Ophthalmoscope Handle (Any Voltage)
- First class laryngoscope blades designed for vet use

*The magnifier on the 3 smaller sizes is amazing and makes the scope useful for; helping to get seeds out, hair stuck in the back of a cats throat, as an otoscope, as an anoscope, as a vaginoscope, as nasal illuminator etc etc

These laryngoscope blades are designed specifically for veterinary use. Based on the popular straight blade Miller style the Opticlar Uni-Fit Blades are available in 00, 0, 1, 2 and 3 sizes.

Code	Description	Blade Size	Magnification	Handle Fit	Price
To Fit Opticlar & WA Handles					
760.020.0000L	Miller Vet	00	No	OPTICLAR & WA	£89.00
760.020.000L	Miller Vet	0	No	OPTICLAR & WA	£89.00
760.020.001L	Miller Vet	1	No	OPTICLAR & WA	£89.00
760.020.002L	Miller Vet	2	No	OPTICLAR & WA	£89.00
760.020.003L	Miller Vet	3	No	OPTICLAR & WA	£89.00
With Magnification					
760.020.0000LM	Miller Vet 3x Mag	00	Yes	OPTICLAR & WA	£105.00
760.020.000LM	Miller Vet 3x Mag	0	Yes	OPTICLAR & WA	£105.00
760.020.001LM	Miller Vet 3x Mag	1	Yes	OPTICLAR & WA	£105.00
To Fit Heine Handles					
760.020.0000LH	Miller Vet	00	No	Heine	£113.00
760.020.000LH	Miller Vet	0	No	Heine	£113.00
760.020.001LH	Miller Vet	1	No	Heine	£113.00
760.020.002LH	Miller Vet	2	No	Heine	£113.00
760.020.003LH	Miller Vet	3	No	Heine	£113.00
With Magnification					
760.020.0000LMH	Miller Vet 3x Mag	00	Yes	Heine	£129.00
760.020.000LMH	Miller Vet 3x Mag	0	Yes	Heine	£129.00
760.020.001LMH	Miller Vet 3x Mag	1	Yes	Heine	£129.00
To Fit Keeler Handles					
760.020.0000LK	Miller Vet	00	No	Keeler	£103.00
760.020.000LK	Miller Vet	0	No	Keeler	£103.00
760.020.001LK	Miller Vet	1	No	Keeler	£103.00
760.020.002LK	Miller Vet	2	No	Keeler	£103.00
760.020.003LK	Miller Vet	3	No	Keeler	£103.00
With Magnification					
760.020.0000LMK	Miller Vet 3x Mag	00	Yes	Keeler	£119.00
760.020.000LMK	Miller Vet 3x Mag	0	Yes	Keeler	£119.00
760.020.001LMK	Miller Vet 3x Mag	1	Yes	Keeler	£119.00

Adaptors only to convert OC/WA blades to Heine or Keeler Fitting:

Code	Description	Price
100-000-099	Heine handle to Blade	£23.99
100-000-097	Keeler handle to Blade	13.99

Any questions please call us on 0800 0683300